## Gebrauchsinformation »Self adhesive composite luting cement«

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## Instructions for use »Self adhesive composite luting cement«

*Self adhesive composite luting cement* is a universal, self adhesive, dual-curing composite luting cement for cementing indirect restorations and endodontic posts. *Self adhesive composite luting cement* does require etching or priming on the tooth.

### Indication

Permanent luting of inlays, onlays, crowns and bridges made of porcelain, oxide porcelain and metal; as well as permanent luting of endodontic posts made of fiber reinforced composite.

### Contraindication

*Self adhesive composite luting cement* is contraindicated, if the preparation does not offer sufficient retention (e.g. Veneers).

*Self adhesive composite luting cement* should not be applied onto the open pulp or dentine close to the pulp.

### Handling of double syringe

The double syringe automatically dispenses and mixes the material, which saves time when applying the luting material directly into the restoration. Remove the transport cap before using a new double syringe and then attach the mixing tip. Following application, the mixing tip acts as a seal and can be left on the double syringe until the next use. The mixing tip should only be replaced by a new tip when the double syringe is used again.

### Recommended use

The following times apply as soon as *Self adhesive composite luting cement* is extruded from the double syringe (working time and curing time depend on the ambient temperature).

- Working time (including mixing time): approx. 1:00 min (room temperature 23°C / 73°F)
- Curing time (including mixing time): approx. 7:00 min (intra-oral)

Caution! After extrusion from the syringe, curing of *Self adhesive composite luting cement* can be accelerated by strong lighting.

### Application

1. Remove the temporary restoration and carefully clean off any residual temporary luting cement from the cavity or prepared tooth. Rinse with water spray.
2. Then check the fit of the restoration.
3. After checking the fit, thoroughly clean the restoration and condition as follows:
  - Etchable porcelain  
Etch preparation with a 5% hydrofluoric etching gel (please follow the manufacturer's instructions). Rinse with water thoroughly and dry with oil-free air. Apply a silane (e.g. *Silane primer*) (please follow manufacturer's instructions).
  - Fiber reinforced endodontic posts  
Clean endodontic post with ethanol, dry with oil-free air and apply silane (e.g. *Silane primer*) (please follow manufacturer's instructions).
  - Other restorations  
Follow manufacturer's instructions for conditioning.
4. Clean cavity / preparation with water spray and dry with oil-free air or pellets. Avoid excessive residual moisture as well as over drying.
5. Apply the *Self adhesive composite luting cement* directly into the restoration using the tip. The working time of *Self adhesive composite luting cement* is approx. 1 min. Seat the restoration fully, applying light pressure. Allow the patient to move into habitual occlusion gradually. A cotton wool roll is a practical aid for this. The restoration can be fully loaded after 7 min.
6. If required, cure material for a minimum of 20 s using a suitable light-curing unit. (The unit should have a minimum light intensity of 400 mW/cm<sup>2</sup>, emit at 450 nm and be regularly checked.)

## Removing excess material

### Self curing

Any excess should be removed during the elastic phase, i.e. approximately 0:30 to 1:00 min after placement, using a micro-brush/ brush/ foam pellet/ dental floss or scaler. It is particularly important to ensure that excess in areas difficult to access (proximal, gingival margins) is removed quickly.

### Self-curing with additional light curing

Light cure the excess cement briefly (1-2 s), which allows it to be easily removed with a scaler. It is particularly important to ensure that excess material in areas difficult to access (proximally, gingival margins) is removed quickly, as it is difficult to remove excess from these areas once the material has fully cured.

### Please note

Do not use *Self adhesive composite luting cement* together with materials that contain eugenol. We therefore recommend eugenol-free temporary cements (e.g. *Temporary luting cement*).

Avoid skin contact with the paste. On accidental contact, wash with soap and water. On contact with the eyes, rinse immediately with plenty of water and consult a doctor.

### Side effects

There are no known systemic side effects. In individual cases a hypersensitivity to components of the material cannot be ruled out.

### Composition

Barium glass in a BIS-GMA-based matrix comprising dental resins, pigments, additives and catalyst. Filler content: 69% by weight = 46% by volume. The variation width of the inorganic filler particles is between 0.02 and 9.3 µm.

### Classification

*Self adhesive composite luting cement* complies with the requirements of ISO 4049:2000.

### Storage

Store cool between 2-8°C / 36-46°F. Do not use after date of expiry.

### Available product sizes

#### Self adhesive composite luting cement

REF 2064 1 Double syringe @ 5 ml

10 Mixing tips

**Keep away from children! For dental use only!**