»VPS Heavy body/VPS Heavy body (fast)«

de

Produktbeschreibung

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Das Produkt entspricht ISO 4823:2021, heavy bodied consistency

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en-US

Product description

VPS Heavy body/VPS Heavy body (fast) is a high consistency precision impression material based on VPS addition-cured silicones and is automatically mixed.

Indication

Dental cases that require the production of dental models.

Contraindication

Do not use the material where there are contact allergies or known allergies to any of the components.

Intended purpose

Impressions using the final impression or double-mixing technique combined with *VPS Light* for the fabrication of crowns, bridges, inlays and onlays.

Limitations of use

Use the material for two-phase impressions only.

Intended user

Dentist, dental assistant staff

Patient target group

People treated in the course of a dental procedure.

Recommended use

Preparing the impression

Have *VPS Light* correction impression material ready for use, and shortly before application, check that the outlets of the safety cartridges are not obstructed. The used mixing tip acts as a seal until the next application is carried out. Do not reuse the original closure cap of the safety cartridge under any circumstances! This can result in unwanted curing due to contamination of the base and catalyst pastes. If retraction fluids have been used, carefully rinse and dry prior to taking the impression.

Corrective impression

Place the impression tray filled with *VPS Heavy body/VPS Heavy body (fast)* into the mouth. Maintain this position for at least the recommended setting time in the mouth for the material. After removal, thoroughly clean and dry the initial impression. Prior to performing the correction impression, cut drainage channels into the initial impression and eliminate any defects.

Fill *VPS Light* correction impression material into the pre-impression, and if necessary, inject it around the preparations. Place the filled final impression back into the mouth, observing the working and setting time of the correction impression material. To ensure that the two impression materials bond securely together, do not wait longer than 30 min between the completion of the initial impression and the preparation of the final impression.

Double-mixing technique

While the assistant loads the impression tray, the dentist can inject *VPS Light* around the prepared teeth. Be sure to fill the impression tray and inject around the prepared teeth as simultaneously as possible, since the increased

temperature of the patient's mouth causes the correction impression material to cure more quickly than the tray material. Place the impression tray into the mouth. Maintain this position for at least the recommended setting time in the mouth for the material.

Disinfection and casting

The impression can be disinfected using commercially available dipping solutions. The use of a 2 % glutaraldehyde solution is recommended. Please ensure that the disinfectant used is indicated for use with addition-cured silicone impression materials.

The exposure time for the impression is approx. 10 to 15 min. Afterwards, rinse the impression under running water for 15 s. Store the impression at normal room temperature (max. 25 °C/77 °F). Wait at least 30 min after removing the impression from the patient's mouth before casting.*

The impression material is compatible with all commercially available plasters. The use of a hard plaster (ISO 6873:2013, Type 3, e.g., Moldano, Bayer) for the cast and an ultra-hard dental plaster for the molds (ISO 6873:2013, Type 4, e.g., Fujirock II, GC) is recommended. The impressions can be copper or silver-plated.

*For more information, see the statement issued by the German Society for Dental and Oral Medicine (DGZMK). Comments of the DGZMK 3/93 V 1.0; www.dgzmk.de ADA recommends following the advice of the Council on Dental Materials, Instruments and Equipment. Disinfection of impressions. JADA 1991;122(3):110; www.ada.org

Storing the impression

- Store the impression at 15 – 25 °C (59 – 77 °F) in a dry place!
- Protect against exposure to direct sunlight!
- The impression is stable and can be cast even after months without loss of accuracy

Time schedule in tabular form

VPS Heavy body	
Working time (from the start of mixing)*	≤ 105 s
Recommended setting time in the mouth	≥ 210 s
VPS Heavy body (fast)	
Working time (from the start of mixing)*	≤ 60 s
Recommended setting time in the mouth	≥ 150 s

* The specified working times are valid for a room temperature of 23 °C (72 °F) and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

The working time describes the period after mixing has begun until the time by which the impression tray must be positioned in the patient's mouth.

Notes for use

- The recommended times can be found in the section headed »Timing«.
- Perforated or non-perforated metal or resin trays can be used with the impression material.
- When selecting the impression tray, ensure that there is enough room between the wall of the impression tray and the dental arch.
- In order to improve the adhesion of the impression material to non-perforated impression trays, apply a thin, even film of *Tray adhesive for VPS*. Allow *Tray adhesive for VPS* to dry for 3 min before loading the impression tray.

»VPS Heavy body/VPS Heavy body (fast)«

Interactions

- Do not use with condensation-curing or polyether impression materials.
- Do not allow the impression to come into contact with organic solvents or liquids that contain solvents as this can cause the material to swell and lead to inaccuracies.
- As retraction solutions may impair the setting time of the material, any retraction solution residues must be removed prior to taking the impression by thoroughly spraying and drying the region.
- Do not use with latex gloves because they inhibit the setting process. The use of PVC or PE-based gloves is recommended.


Cartridge handling

Insert the cartridge into the mixing device. Please refer to the operating instructions of the respective device manufacturer.


Attaching/changing the mixing tip

Open the locking clip by pulling it along the track to its terminal position and pull off the original closure cap or the old mixing tip. (Fig 1-2)

1


2


Then insert a new mixing tip and press both openings on the back of the mixing tip firmly onto the corresponding openings of the cartridge. (Fig 3)

3

Lock the attached mixing tip by sliding the locking clip along the track until it encloses the mixing tip and reaches its final position. (Fig 4)

4

Note: Always replace the mixing tip before each mixing procedure! The mixers are unable to press out impression material through material hardened in a used mixing tip.

Note: Make sure that the mixing tip sits firmly on the cartridge! Otherwise, the impression material will leak from the cartridge and mixing tip during extrusion, which may cause contamination.

Note: Prior to using a new cartridge for the first time, discard approx. the first 3 cm of the ejected material. After use, leave the used mixing tip on the cartridge as a seal.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

»VPS Heavy body/VPS Heavy body (fast)«

Product performance characteristics

The product complies with ISO 4823:2021, heavy bodied consistency

VPS Heavy body

Mixing ratio of base/catalyst 5:1

VPS Heavy body (fast)

Mixing ratio of base/catalyst 5:1

Storage/disposal

- Store in a dry place at 15 to 25 °C (59 to 77 °F)!
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

REF 9048MDR: SiO₂, vinyl silicones, hydrogen silicone, additives, pigments, catalysts.

REF 9049MDR: SiO₂, vinyl silicones, hydrogen silicone, additives, pigments, catalysts.

Packaging

VPS Heavy body	
REF 9048MDR	1 Cartridge @ 380 ml, 10 Mixing tips
VPS Heavy body (fast)	
REF 9049MDR	1 Cartridge @ 380 ml, 10 Mixing tips