## »Adhesive for composite«

de

### Produktbeschreibung

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## »Adhesive for composite«

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EBPADMA, mmHE (ca. 99 %), Additive, Katalysatoren

Dentalharze (Bis-GMA, TEDMA, UDMA, EBPADMA; ca. 99 %), Additive

### Handelsformen

REF 2051MDR	1 Flasche à 5 ml
	10 Pinsel
	1 Mischpalette

## »Adhesive for composite«

en-us

### Product description

*Adhesive for composite* is a light-curing bonding agent for the intra-oral repair of damaged composite restorations. It bonds newly-applied composite restoration materials to already-cured older composites. Use of this intra-orally eliminates the need for fully removing the restoration.

### Indications

Loss of tooth substance due to caries, trauma, tooth attrition/abrasion or development disorders

### Contraindications

Do not apply material directly onto exposed dentin surfaces or onto the exposed pulpa. Do not use the material where there are contact allergies or known allergies to any of the components.

### Intended use

Intra-orally repairing damaged composites, correcting the shades of esthetically unsatisfactory composite restorations.

### Limitations of use

The use of *Adhesive for composite* is contraindicated if a dry operating field cannot be ensured during use.

### Intended users

Dentist

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Recommended use

1. Roughen up the restoration on the surface to be repaired with a coarse-grit grinding disc or a diamond bur (Ra approx. 150 µm) and round off the edges of the fracture. It is advisable to use a rubber dam.
2. Etch exposed enamel surfaces with a 37% phosphoric acid *etching gel* (e.g. *Etching gel*) before applying *Adhesive for composite*. Observe the manufacturer's instructions for the *etching gel* used. Exposed dentin surfaces must be treated with an adhesive system (e.g. *Self etching adhesive*). In doing so, it is essential to observe the instructions from the manufacturer of the adhesive system.
3. Use a brush to apply a drop of *Adhesive for composite* onto the surface to be repaired and work in for at least 30 seconds. When repairing chemically cured composites, this time increases to approx. 60 seconds.

**Important:** These times must be complied with under all circumstances so that the adhesive can bond properly to the restorative.

4. Gently blow the material with oil- and water-free air so that a slightly damp layer remains on the surface.
5. Then light-cure the *Adhesive for composite* for at least 20 seconds. A light-curing time of 30 seconds is recommended when using chemically cured composites.

6. Use the restorative (e.g. *Nano hybrid composite*) in accordance with the instructions of the manufacturer in question.

### Notes for use

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- Do not apply material directly onto exposed dentin surfaces.
- Light units should have a light intensity of 600 mW/cm<sup>2</sup>, emit at 450 nm, and their performance should be checked regularly. The light-curing unit should be positioned as closely to the material as possible.
- The oxygen inhibition layer that develops during use facilitates bonding with the materials applied in subsequent layers, and therefore must not be removed.
- Contamination of the bottle's contents with a contaminated brush and premature polymerization due to ambient light should be avoided.

### Interactions

*Adhesive for composite* should not be used in conjunction with materials that contain eugenol, as these can inhibit polymerization.

### Residual risks / side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

### Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Performance characteristics of the device

SBS bonding properties on cured composite materials for repair are ≥ 10 MPa

### Storage / disposal

- Dry at room temperature (15-25 °C/59-77 °F) and store in a location that is protected against light.
- Do not use after the expiration date.
- Disposal must comply with national regulations.

### Composition

EBPADMA, mmHE (approx. 99 %), additives, catalysts

Dental resins (bis-GMA, TEDMA, UDMA, EBPADMA; approx. 99 %), additives

**»Adhesive for composite«****Packaging**

REF 2051MDR	1 Bottle @ 5 ml
	10 Brushes
	1 Mixing pad