## »Composite core build-up z dual«

de

### Produktbeschreibung

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en-US

### Product description

*Composite core build-up z dual* is a dual-cure, radiopaque core build-up composite that is automatically mixed and directly applied. Using the Intraoral-Tips provided, even difficult regions are easily accessed, such as in root post cementation.

### Indications

Severe loss of dental hard tissue

### Contraindications

Do not use the material where there are contact allergies or known allergies to any of the components.

### Intended use

- All types of core build-ups and build-up fillings
- Cementation of root posts

### Limitations of use

Do not use directly on the exposed pulp

### Intended users

Dentist

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Recommended use

#### Core build-up

1. To facilitate application of *Composite core build-up z dual*, place a matrix if necessary around the core that is to be prepared.
2. Apply *Composite core build-up z dual* directly inside the preparation. The material can be formed using commercially available instruments, e.g. a Heidemann spatula. The working time for *Composite core build-up z dual* is 01:30 min. Curing can be accelerated by exposure to a suitable dental curing light. The maximum layer thickness should be 4 mm if an exposure time of at least 40 seconds is used; with an exposure time of at least 20 seconds, the maximum layer thickness should be 2 mm. With a chemical self-cure only, the restoration must be allowed to cure for at least 06:00 min before it can support a full load.
3. Remove the matrix after *Composite core build-up z dual* has cured. If a pre-fabricated core build-up form is used, *Composite core build-up z dual* can be applied directly into the core build-up form.
4. Then place the core build-up form onto the prepared stump and keep in position until *Composite core build-up z dual* has cured.

#### Cementation of root posts

1. Prepare the root post in accordance with the instructions for use provided by the manufacturer.
2. Following preparation of the root canal, etch the inside of the canal as per the instructions of the manufacturer of the bonding agent used.

3. Work in the bonding agent. Use a paper tip to remove any excess prior to light curing.
4. Apply *Composite core build-up z dual* inside the prepared canal. The root post can be coated with *Composite core build-up z dual* before it is placed into the filled canal.
5. Once the root post has been positioned, light cure the material using a suitable lamp for 20 seconds in order to secure the root post. As *Composite core build-up z dual* is dual-curing, areas that cannot be reached by light will self-cure chemically within 06.00 min. *Composite core build-up z dual* can also be used for the core build-up of the coronal portion of the tooth (see the section on "Core build-up"). A matrix should be used here due to the low viscosity of the material.

### Fabricating a temporary restoration

If a resin-based temporary restoration is to be placed onto the prepared stump, use a suitable separation medium (e.g. Vaseline) for isolation.

### Timing

Application inside the preparation	00:00 min
End of working time (chemical curing only)	01:30 min
Curing with max. 2 mm layer thickness	Light exposure for 20 seconds
Curing with max. 4 mm layer thickness	Light exposure for 40 seconds
<i>Composite core build-up z dual</i> is chemically cured	06:00 min

### Notes for use

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- **Cavity preparation**

**Note:** It is advisable to use a rubber dam.

1. Before applying *Composite core build-up z dual*, remove caries and existing fillings and provide for appropriate pulp/dentine protection.
2. Irrigate the cavity carefully and air dry. Avoid both excessive residual moisture and overdrying. For adhesive anchorage of *Composite core build-up z dual* to the hard tissue, a dual-cure bonding agent should be used (e.g. *Total-etch adhesive dual*). Observe the instructions for use provided by the manufacturer.

- Light-curing units should emit at 450 nm and their light output should be checked regularly. The luminous intensity should be at least 600 mW/cm<sup>2</sup>. Place the light-curing unit as close as possible to the material.
- The material should be at room temperature before use. Cooler material cures more slowly and has a higher viscosity.

### Interactions

- The use of *Composite core build-up z dual* with single-bottle bonding systems is contraindicated.

## »Composite core build-up z dual«

- Do not use *Composite core build-up z dual* together with materials that contain eugenol!

### Using the double-barreled syringe

The double-barreled syringe enables automatic dispensing and mixing of the material, which allows a time-saving direct application. When using a new double-barreled syringe, remove the transport cap first. The mixing tip can then be fitted. After use and until the next treatment, leave the mixing tip on the double-barreled syringe as a seal. The used mixing tip is replaced with a new tip when the device is used for the next treatment.

---

**Note:** When using a double-barreled syringe for the first time, discard a small amount of the initially ejected material (e.g. a pea-sized amount).

---

### Residual risks / side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

### Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Performance characteristics of the device

The product complies with ISO 4049:2019, type 2, class 3.

### Storage / disposal

- Store in a dry place at 2- 25°C/35.6- 77°F!
- After application, leave the used mixing tip on the cartridge or syringe as a cap!
- Do not use after the expiration date.
- Disposal must comply with national regulations.

### Composition

Dental glass, silica (sil), bis-GMA, UDMA, HDDMA, DDDMA, mmHE, additives.  
Inorganic filling materials: approx. 43 vol% (0.02 – 2.5 µm).

### Packaging

REF 2100MDR	1 Double-barreled syringe (A3) @ 5 ml
	10 Mixing-Tips
	10 Intraoral-Tips