## Gebrauchsinformation »Nano hybrid flowable composite«

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### Instructions for use »Nano hybrid flowable composite«

*Nano hybrid flowable composite* is a light cure, radiopaque filling material for the restorative therapy.

#### Indication

As initial layer under Class I and II restorations, Anterior restorations (Classes III, IV), Class V restorations (cervical caries, root erosion, wedge-shaped defects), Small restorations of all types, Blocking out of undercuts, Repair of composite and ceramic veneers, Repair of temporaries

#### Recommended use

1. Select the shade of the filling composite
2. Provide appropriate isolation, best with a rubber dam
3. Prepare the cavity according to the requirements of the adhesive technique
4. Clean the cavity with water spray
5. Air dry the cavity
6. Apply a pulp protector (calcium hydroxide) if needed; only cover areas close to the pulp and subsequently apply a pressure-resistant cement
7. Apply matrix / interdental wedge
8. Condition and apply the bonding agent according to the instructions for use of the product in use
9. Apply *Nano hybrid flowable composite* in layers of max. 2 mm
10. Polymerize each layer for 20 s using a light intensity of at least 500 mW/cm<sup>2</sup>. With a light intensity of 1100 mW/cm<sup>2</sup>, the individual layers can be cured in 10 s. When using a metal matrix (e.g. OptraMatrix), additionally polymerize from the buccal or lingual/palatal after removing the matrix
11. Finish the restoration with suitable finishers or fine diamonds
12. Check the occlusion
13. Polish with silicone polishers as well as with polishing discs and polishing strips

#### Please note

- Light curing units should have 450 nm output and should be checked regularly. The light intensity should be at least 500mW/cm<sup>2</sup>. Place the light as close as possible to the restorative material.
- Do not remove the oxygen-inhibited layer as it is needed for bonding to the next layer.
- Follow the instructions for use of any materials, which are used with *Nano hybrid flowable composite*.

#### Interactions

Avoid materials containing eugenol, moisture or oily air; they can inhibit the polymerization at the contact area. The use of cationic mouthwashes, plaque indicators or chlorohexidine may lead to discolorations.

#### Contraindications

Application of *Nano hybrid flowable composite* is contraindicated, when a dry working area or the recommended application technique is not possible, for direct application on the opened pulp or when allergies to any components of *Nano hybrid flowable composite* are known.

#### Side effects

Side effects are not known to date. In singular cases, hypersensitivity to components of the material cannot be excluded. Avoid contact of paste with skin. If accidental contact occurs, wash immediately with soap and water. In case of contact with eyes, rinse immediately with plenty of water and consult a physician.

#### Composition

Barium glass in a Bis-GMA-based matrix of dental resins. pigments, additives and catalyst. Filler ratio: 65 percent by weight = 38 percent by volume. The variation of anorganic filler particle is between 0,02-3 µm.

### Classification

*Nano hybrid flowable composite* fulfills the requirements of ISO 4049:2009.

The stated value of the depth of cure according ISO 4049:2009 of *Nano hybrid flowable composite* is 2 mm. It is recommended, in any case, to apply the material in layers of 2 mm and to light cure each layer.

#### Storage

Store at normal room temperature (15-25 °C / 59-77 °F). Please close syringe immediately after use to prevent premature setting owing to ambient light. Do not use after date of expiry.

#### Available product sizes

##### Nano hybrid flowable composite

REF 2060 2 Syringes (A2) @ 1.5 g

10 Luer lock tips

**Keep away from children! For dental use only!**