## »Self etching flowable composite«

de

### Produktbeschreibung

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en-US

### Product description

*Self etching flowable composite* is a self-etching, self-adhesive, radiopaque, flowable composite that is cured with light. It is immediately ready for use because the preparatory steps of etching, priming and bonding are contained in it.

### Indications

- Loss of tooth substance due to caries, trauma, tooth attrition/abrasion or development disorders
- Molars and premolars with pronounced deep fissures and increased caries risk: Fissure sealing

### Contraindications

- Do not apply the material to the exposed pulp.
- Do not use the material where there are contact allergies or known allergies to any of the components.

### Intended use

- Small restorations of class I and III
- Restorations of class V
- Underfilling material for class I and II
- (Extended) fissure sealing
- Repairs of composite restorations
- Modifications to temporaries and long-term temporaries
- Blocking out and filling of undercuts in cavities
- Small occlusal primary tooth cavities

### Limitations of use

- Do not use the material if a sufficiently dry working area or the recommended application technique is not possible.
- Do not use the material on etched dentin

### Intended users

Dentist

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Recommended use

#### Small restorations of class I, III and V and underfilling material for class I and II

1. Determine tooth color and select the appropriate shade.
2. Prepare the cavity in accordance with the general rules for adhesive technique. Bevel the enamel edges if necessary. Protect pulp in deep cavities with appropriate preparations.

**Please observe the following:** It is advisable to use a rubber dam.

3. Clean tooth, remove all cleaning residue with water spray.

**Please observe the following:** Cleaning with alcohol may cause overdrying of the dentin.

4. Dry the tooth with water-free and oil-free air. Avoid both excessive residual moisture and overdrying. A moist layer must remain on the surface of the tooth.
5. Apply the material to the cavity surface with the aid of the Luer-Lock-Tips by pressing the syringe and massage a thin layer ( $\approx 0.5$  mm) into the entire surface of the cavity wall for 25 seconds using the brush.
6. Remove any excess, if necessary.
7. Polymerize for 20 seconds with a light intensity of  $1000 \text{ mW/cm}^2$  and for 40 seconds with a light intensity of  $600 \text{ mW/cm}^2$ .
8. Layer the material as follows:
  - Small restorations of class I, III and V: Build up layers with a maximum thickness of 2 mm with the material and light-cure each layer for 20 seconds at a light intensity of  $1000 \text{ mW/cm}^2$  and 40 seconds at a light intensity of  $600 \text{ mW/cm}^2$ .
  - Underfilling material for class I and II: Build up a layer with a maximum thickness of 1 mm with the material and light-cure for 20 seconds at a light intensity of  $1000 \text{ mW/cm}^2$  and 40 seconds at a light intensity of  $600 \text{ mW/cm}^2$ .

**Please observe the following:** It is not necessary to work in the individual layers once more.

9. Remove excess, contour the filling surface and finish using suitable instruments.
10. Use a polishing paste to polish the filling surface and neighboring surfaces. Use a polishing cup or polishing disc for polishing work.

#### Fissure sealing and extended fissure sealing

1. For extended fissure sealing, widen particularly narrow and vial-shaped fissures.
2. Clean enamel with fluoride-free paste, remove all cleaning residue with water spray and dry the tooth with water-free and oil-free air. In doing so, leave the enamel surface moist.
3. Determine tooth color and select the appropriate shade.

**Please observe the following:** It is advisable to additionally etch the enamel. To do this, etch the enamel for 15 – 30 seconds with 37% phosphoric acid gel, rinse for 15 seconds with water and dry with oil-free and water-free air. In doing so, leave the enamel surface moist.

4. Apply the material to the surface of the fissure or cavity wall with the aid of the Luer-Lock-Tips by pressing the syringe and massage a thin layer ( $\approx 0.5$  mm) into the entire surface for 25 s using the brush.
5. Remove any excess, if necessary.
6. Light-cure for 20 or 40 seconds depending on the light intensity of the lamp.
7. Where necessary, use the material to build up layers with a maximum 2 mm layer thickness and light-cure for 20 or 40 seconds.

## »Self etching flowable composite«

**Please observe the following:** It is not necessary to work in the individual layers once more.

8. Check the occlusion and correct if necessary.

### Notes for use

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- Light-curing units should have a luminous intensity of at least 600 mW/cm<sup>2</sup>, emit at 450 nm, and their performance should be checked regularly. Observe the manufacturer's instructions for the light-curing unit used.
- The light-curing unit should be positioned as closely to the material as possible.
- It is advisable to use a rubber dam.
- Do not remove the oxygen inhibited layer because this is required for the bond with the next layer.
- Observe the instructions for use of the materials that are used with the *self etching flowable composite*.

### Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

### Residual risks / side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

### Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- The brush holder is reusable after disinfection, although it should be replaced if visible signs of wear, e.g. scratches, appear. Disinfection can be carried out using a commercially available wiping or dipping disinfectant. It is advisable to only use disinfecting solutions listed by the RKI (Robert Koch Institute).
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Performance characteristics of the device

The mechanical properties of the product comply with ISO 4049: 2019, Type 1, Class 2, Group 1.

Polymerization depth  $\geq 2$  mm with a 40 second light-curing time when the light intensity of the lamp is  $\geq 600$  mW/cm<sup>2</sup>. Polymerization depth  $\geq 2$  mm with a 20 second light-curing time when the light intensity of the lamp is  $\geq 1000$  mW/cm<sup>2</sup>

### Storage / disposal

- Store in a dry place at 2- 25°C/35.6- 77°F!
- Reseal syringes after use.
- Do not use after the expiration date

### Composition

Dental glass, bis-GMA, EBPADMA, alip. EDA, UDMA, MDP, HEMA, TEDMA, additives.

Inorganic filling materials: approx. 43 vol% (0.7 – -2.3  $\mu$ m).

### Packaging

REF 2109MDR	2 Syringes @ 2 g Composite (A2) 1 Brush holder 20 Brushes 20 Luer-lock tips
-------------	--