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Instructions for use »Etching gel«

Etching gel is a 37% phosphoric etching gel with optimal flow properties for etching tooth structure, e.g. for fissure sealing, in the adhesive technique with direct filling treatment using composites or when luting laboratory-fabricated restorations or bonding orthodontic retentive elements.

Indication

Enamel and dentin etching

Recommended use

Preparation

The areas to be etched should be clearly isolated.

1. The use of a rubber dam is recommended! Avoid any contact with the soft tissue (skin, oral mucosa, eyes etc.).
2. Clean the tooth using a polishing cup with fluoride-free and grease-free polishing paste.
3. Carefully remove the cleaning residue using water and then clean and dry the tooth surface with oil-free and water-free compressed air.
4. Prepare the cavity accordingly; bevel the enamel margins to produce durable filling margins. Rinse the cavity thoroughly using water and dry with oil-free and water-free compressed air.

Protect the pulp in deep cavities (e.g. cover with a thin layer of a suitable cavity liner material). The cavity should be dry and clean.

Etching procedure

1. First, apply *Etching gel* to the enamel margins of the cavity and then apply to the dentin if required.

Recommended reaction time

20 - 60 s on enamel

≤ 15 s on dentin

50 - 80 s on teeth with high fluoride content

2. Then rinse for a minimum of 20 s with water.
3. Dry with oil-free and water-free compressed air.

The etched surface should not be touched or contaminated with saliva before the next stage of treatment. If there is contamination after etching and drying, however, the surface should be briefly etched again (approx. 10 s). The etched enamel should have a chalky white appearance! Adhesives (e.g. *Total-etch adhesive*) etc. can be applied according to the manufacturer's instructions after etching.

Please note

- *Etching gel* causes burns.
- Avoid contact with the skin, eyes and soft tissue. If there is contact, rinse with plenty of water. Consult a doctor if necessary.

Interactions

Cavity lining material may be etched by the etching gel with a longer reaction time.

Side effects

There are no known systemic side effects. In individual cases a hypersensitivity to components of the material cannot be ruled out.

Composition

37% phosphoric acid, water, aerosil.

Storage

Store at normal room temperature (15-25 °C / 59-77 °F). Do not use after date of expiry.

Available product sizes

Etching gel

REF 2052 1 Syringe @ 2 ml

Keep away from children! For dental use only!