Gebrauchsinformation »Composite core build-up z dual«

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Instructions for use "*Composite core build-up z dual*"

Composite core build-up z dual is a dual-cure, radiopaque core build-up composite that is automatically mixed and directly applied. Using the intra-oral tips provided, even difficult regions are easily accessed, such as in root post cementation.

Indications

- All types of core build-up and build-up fillings
- Cementation of root posts

Using the double-barreled syringe

The double-barreled syringe allows automatic dosage and mixing of the material.

When using a new double-barreled syringe, remove the transport cap. Attach the mixing tip and an intra-oral tip if necessary.

When using for the first time, discard the material that is initially expressed from the mixing tip (approx. a pea-sized amount). A uniform mixing result is then guaranteed.

After use, the mixing tip acts as a cap and can be left on the double-barreled syringe until the next time it is used. Only when the double-barreled syringe is used again must the mixing tip be replaced by a new one.

Cavity preparation

Hint: The use of a rubber dam is recommended.

1. Before applying *Composite core build-up z dual*, remove caries and existing fillings and provide for appropriate pulp/dentine protection.
2. Irrigate the cavity carefully and air dry. Avoid excessive residual moisture and over-drying. For adhesive anchorage of *Composite core build-up z dual* to the hard tissue, a dual-cure bonding agent should be used (e.g. *Total-etch adhesive dual*). Observe the instructions for use provided by the manufacturer.

Recommended use

Core build-up

1. To facilitate application of *Composite core build-up z dual*, place a matrix if necessary around the core that is to be prepared.
2. *Composite core build-up z dual* can then be applied directly inside the preparation. The material can be formed using commercially available instruments, e.g. a Heidemann spatula. The working time for *Composite core build-up z dual* is 1:30 min. Curing can be accelerated through application of light using a suitable lamp. With an application time of at least 40 seconds, the maximum layer thickness is 4 mm; with an application time of at least 20 seconds, the maximum layer thickness is 2 mm. With chemical curing only, the restoration must be cured for at least 6 minutes before it can support a full load.
3. The matrix may not be removed until *Composite core build-up z dual* has cured. If a pre-fabricated core build-up form is used, *Composite core build-up z dual* can be applied directly into the core build-up form.
4. Then place the core build-up form onto the prepared core and keep in position until *Composite core build-up z dual* has cured.

Cementation of root posts

1. Prepare the root post in accordance with the instructions for use provided by the manufacturer.
2. Following preparation of the root canal, etch inside the canal as per the instructions of the manufacturer of the bonding agent used.
3. Work in the bonding agent. Use a paper tip to remove any excess prior to the application of light.
4. Apply *Composite core build-up z dual* inside the prepared canal. Before the root post is set, it can be coated in addition to the filled canal with *Composite core build-up z dual*.
5. Once the root post has been positioned, light cure the material using a suitable lamp for 20 seconds in order to secure the root post. As *Composite core build-up z dual* is dual-curing, areas that cannot be reached by light cure chemically within 6 minutes. *Composite core build-up z dual* can also be used for core build-up of the coronal portion of the tooth (see

the section on "Core build-up"). A matrix should be used here due to the low viscosity of the material.

Fabricating a temporary restoration

If a plastic temporary restoration is to be applied to the prepared core, use a suitable separation medium (e.g. *Vaseline*) for isolation.

Working times

0 minutes	Application inside the preparation.
1:30 minutes	End of working time (chemical curing only).
20 seconds: application of light	Curing with max. 2 mm layer thickness.
40 seconds: application of light	Curing with max. 4 mm layer thickness.
6 minutes	<i>Composite core build-up z dual</i> is chemically cured.

Please note

- Light-curing units should emit at 450 nm and their light output should be checked regularly. The luminous intensity should be at least 400 mW/cm². Place the light as close as possible to the material.
- The use of *Composite core build-up z dual* with single-bottle bonding systems is contraindicated.
- Do not use *Composite core build-up z dual* together with materials that contain eugenol.
- The material must be at room temperature for use. Cool material cures more slowly and has a higher viscosity.
- Do not use *Composite core build-up z dual* in cases of known allergies to any of the components, or in the rare event of contact allergies.
- Avoid skin contact with the paste. In case of inadvertent contact, rinse with water and soap. In case of contact with the eyes, rinse immediately with copious amounts of water and consult a physician.

Side effects

No systemic side effects are known to date. Hypersensitivity to components of the material cannot be excluded in individual cases.

Composition

Barium glass, zirconium dioxide in a Bis-GMA-based resin matrix, pyrogenic silicic acid, initiators (including benzoyl peroxide), co-initiators, stabilizers and additives. Filler content: 64 wt % = 43 vol. % (0.02-2.4 µm).

Classification

Composite core build-up z dual in the double-barreled syringe meets the requirements of ISO 4049:2009, type 2, class 3.

Storage

Store at 2-25 °C (36-77 °F). It is recommended that material that has not been opened be stored in a fridge. Do not use after the expiration date.

Packaging

REF 2100 1 double-barreled syringe (A3) with 5 ml

10 mixing tips

10 intra-oral tips

Keep out of reach of children! For dental use only!